LEVANTAR Floor Lift, Power Base 500 lbs



HCPCS: E0635


FIGURE A


FIGURE B


FIGURE C


FIGURE D


FIGURE E



FIGURE F


FIGURE G


FIGURE H

- Weight capacity of 500 lbs
- Base option: Power
- Elegant, contemporary design with soft lines but still heavy duty
- Features you expect in a Premium lift
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height
- offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS	Description
FLP500	Levantar Floor lift, power base, 500 lbs

SPECIFICATIONS	Description
Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Battery Warranty	1 yr
Boom Height	33" - 76"
Warranty on Battery	1 yr

ACCESSORIES	Description
FLSCALE	Patient Lift Scale 660lbs capacity, EA
13012	Patient Sling,Canvas,1/bx, BX
13019	Chains for Patient Lft,1pr/bx, BX
13025	Sling,Polystr.for Lift 1/ea, BX
13026	Sling,Polystr.w/Cmde Open1/ea, BX