

Pinniped Pediatric Commode

HCPCS: E1399-1,E1399-1



FIGURE A



FIGURE B

- All-in-one welded steel frame available in red and blue
- Plastic seat and contoured back
- Durable snap-on toilet seat
- Height adjustable seat allows child to firmly place their feet on the floor for comfort, safety and security
- Adjustable hook & loop chest strap for additional support
- Molded plastic armrests provide additional comfort
- Rubber non-skid tips ensure that the legs will not mark the floor
- Complete with 7.5 qt removable commode bucket with metal handle and cover

PRODUCTS	Description
PC 1000 BL	Pinniped Pediatric Commode, Blue, 1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	250 lbs
Carton Shipping Weight	13 lbs
Overall Width	17"
Weight	9.5 lbs