Bariatric Foam Mattresses

HCPCS: E0184-1,E0184-1






FIGURE A

• Polyurethane Stretch cover is fluid-resistant and easy to clean for extended life and service
• Fire Retardant: Mattress meets the requirements of 16 CFR Part 1632 Standard for the Flammability of Mattresses and Mattress Pads and 16 CFR Part 1633 (Federal Flammability (open flame) standard for mattress sets) when used without a foundation
• High-Density foam mattress for the bariatric patient
• Solid foam inner core provides comfort, support and durability in a one-piece construction
• Latex free

PRODUCTS	Description
15301	Bariatric Mattress, 42in,1/ea
15301-84	Bariatric Foam Matt,84"x42"x6.5",1/ea

SPECIFICATIONS	Description
Warranty	1 Year
Product Dimensions	15301: 80(L) x 42(W) x 6.5(H)|15301-84: 84(L) x 42(W) x 6.5(H)|15310-84: 84(L) x 48(W) x 6.5(H)
Weight Capacity	15301: 585 lbs|15301-84: 585 lbs|15310-84: 665 lbs
FireRating	CFR 16 Part 1633 Federal Fire Code
Side Rails	No
Cover Material	Polyurethane Stretch
Limited Warranty	1 Year
Mattress Layers of Foam	1

ACCESSORIES	Description
15301-84CVR	Cover Bari Foam Matt,84"x42"x6.5",1/ea, EA
15310-84CVR	Cover Bari Foam Matt,84"x48"x6.5",1/ea, EA
15312-84CVR	Cover Bari Foam Matt,84"x54"x6.5",1/ea, EA
15301C	Bariatric Mattress Cover, 42", 1/ea, EA
15310C	Mattress Cover 48" 1/cv, CS