

View SpO2 Deluxe Pulse Oximeter

HCPCS:



FIGURE A



FIGURE B

- Dual color high resolution east-to-read OLED display, auto-rotation multi-display screen
- Display SpO and Pulse Rate, Pulse bar and Plethysmogram
- Now with Perfusion Index and programmable alarm functionality
- Low power composition, Auto-off function
- Approximately 30 hours of normal operation
- Accommodates widest range of finger sizes from pediatric to adult
- Includes lanyard, carry case, and two AAA batteries
- Perfect for sports enthusiasts to measure oxygen saturation and pulse rate

PRODUCTS	Description
MQ3200	View Sp02 Fingertip Pulse Oximeter 1/ea

SPECIFICATIONS	Description
Warranty	Limited 2 Year
Battery Duration	20 to 30 hours
Battery Type	2 AAA alkaline or rechargeable
Dimensions	2.2(L) x 1.2(W) x 1.3(H)
Weight	2.0 oz