PRISM
ORION BACK





The Orion Back ™ was designed for users who require mild to moderate postural support. Extras,such as the positioning options, the ease of attachment and the unique method of adjusting for a wide range of angles, makes the Orion Back appealing for a wide variety of users.

PRODUCT HIGHLIGHTS

- Lightweight, sturdy 5052 grade aluminum shell
- A full 5" range of forward and rear adjustable
- Fully Adjustable for height, depth and angle of recline
- Attractive cover design
- Adjustable Angle Hardware
- Replaceable foam
- Hardware fits 3/4". 7/8" and 1" tubing and can be installed within a few minutes with the right tools
- **Product Weight:** 4.75 lbs, (18"W x 16"H)
  User Weight: 250 lbs
  Warranty: 3 years limited on back, 90 days on cover, and 1 year on hardware


Prism Series Backs STYLE: **Orion Back**

Stock Sizes – Same Day Shipment

Part Number	Size W x H	Part Number	Size W x H
OB1616	16 " x 16"	OB1820	18" x 20"
OB1620	16" x 20"	OB2016	20" x 16"
OB1816	18" x 16"	OB2020	20" x 20"


Non-Stock Sizes – Shipment within 3 business days

Part Number	Size W x H	Part Number	Size W x H
OB1414	14" x 14"	OB1818	18" x 18"
OB1618	16" x 18"	OB1916	19" x 16"
OB1716	17" x 16"	OB1920	19" x 20"
OB1720	17" x 20"		





Lateral Dimensions - Orion Back		
Heights 11-26"		
Chair Width	DIM B (")	DIM A (")
14"	11.85	5
16"	13.85	5
18"	15.85	5
20"	17.85	5
22"	19.85	5



Cover and Foam Layers

FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6

1-888-737-4011

www.futuremobility.ca

USA

1 Buffalo River Place, Buffalo,
NY, 14210

716-783-9236

www.futuremobility.com