

AquaSense Knurled Chrome Grab Bar with Rotating Flange

HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

Knurled Chrome Grab Bars with Rotating Flange by AquaSense® are perfect where extra support is required in the bathroom. They help prevent bathroom accidents such as slipping, sliding or falling. This knurled chrome grab rail has rotating flanges for easy installation on the side of the bathtub wall, the shower wall or around the toilet and sink area walls. The rotating flanges of the grab bar facilitate the installation; the grab bar can be mounted vertically, horizontally or diagonally. AquaSense® Grab Bars are rust resistant and come complete with instructions and stainless steel hardware.

- Rotating flanges allow the grab bar to be mounted vertically, horizontally or diagonally
- Rust-resistant chrome finish
- Knurled grip provides added security
- Supports up to 250 lbs
- Stainless Steel Hardware Included
- Clearance from wall: 2 1/2"

PRODUCTS	Description
785-130	AqS GRAB BAR 12" ROTATE, 1EA, 12CV, Re
785-132	AqS GRAB BAR 16" ROTATE, 1EA, 12CV, Re
785-134	AqS GRAB BAR 18" ROTATE, 1EA, 12CV, Re
785-136	AqS GRAB BAR 24" ROTATE, 1EA, 12CV, Re
785-138	AqS GRAB BAR 32" ROTATE, 1EA, 12CV, Re