

Padded Transfer Bench

HCPCS: E0247


FIGURE A



FIGURE B


FIGURE C

- Comfortable cushioned seat and backrest
- Tool-free assembly of back, legs and arm
- New pinch-free cover allows for push pins to be depressed without pinching fingers
- A-frame construction provides stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs

PRODUCTS	Description
12005KD-1	Transfer Bench, Padded, KD, 1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	24" (W) x 16" (D)
Weight Capacity	400 lb
Carton Shipping Weight	14 lb|26.5 lb
Weight	11.75 lb

ACCESSORIES	Description
12011SL	SuctionCups,Lg.12011,1pr/bx, BX
12011SS	SuctionCups,Sm.12011,1pr/bx, BX