

Plastic Tub Rail

HCPCS: E0246


FIGURE A



FIGURE B

- Attractive, ergonomically-designed tub rail provides a comfortable, safe grip
- Handle is perpendicular to tub
- 16.5" in height
- Tool-free, easy to install
- Width adjusts from 3" - 6.25"
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
12044	Tub Rail,16 1/2"Plastic,1/bx, RTL

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lb
Height	16.5"
Carton Shipping Weight	7 lb
Dimensions	15" (L) x 4"(W) x 12.5" (H)
Weight	5 lb