TRAVEL BUGGY

See *Your* World Again





CITY 2 PLUS

400 lbs
Weight capacity





All-terrain power

115°
Reclining back



KEY FEATURES

POWER THAT GOES FURTHER

High-performance motors deliver smoother rides, improved speed, and extended battery life.



CAREGIVER REMOTE PAIRING



ALL-DAY COMFORT

20.5" tall padded backrest and 3" thick seat cushion provide ultimate support.



3 LIGHTWEIGHT PIECES

Quick-release wheels reduce the weight of the heaviest component to 38 lbs for easy lifting.



DUAL BATTERY RANGE

Dual batteries deliver up to 25 km / 15.5 mi of range—no swapping needed.



BUILT-IN STORAGE



WIDE SEATING ROOM

Install the included armrest spacers to expand the width between armrests to 23".

EXTENDABLE FOOTREST

Stretch out your legs with the 3-position extendable footrest.



SPECIFICATIONS

SKU (BLACK)	22MNA0000109
SKU (SILVER)	22MNA0000107
NET WEIGHT	60 lbs
GROSS WEIGHT	65 lbs
WEIGHT CAPACITY	400 lbs
DRIVING RANGE	25 km / 15.5 mi
INCLINE CAPABILITY	12°
TURNING RADIUS	26"
MAXIMUM SPEED	6 kph / 3.7 mph
MOTOR POWER	2 x 250 W Brushless
UNFOLDED DIMENSIONS	44"L x 23"W x 40"H
FOLDED DIMENSIONS	32"L x 23"W x 13"H
FRONT / REAR WHEEL DIAMETER	7" / 12"

SEATING DIMENSIONS

SEAT SIZE	17.5"W x 19"D x 3"H
BACKREST HEIGHT	20.5"
BETWEEN ARMRESTS	19.5" - 23"
SEAT TO FLOOR HEIGHT	19.5"

BATTERY

BATTERY TYPE	Li-Ion 24 V / 6 Ah
BATTERIES INCLUDED	Two (2)

POPULAR ADD-ONS



HEAVY-DUTY TRAVEL BAG
SKU: 99ASG0400100



ADJUSTABLE HEADREST
SKU: 99ACM0200100



CITY 2 PLUS - CAREGIVER REMOTE
SKU: 22AET0650100



ORTHOPEDIC ELEVATING LEG RESTS
SKU: 22ACM0400100



LIGHTWEIGHT EXTENDABLE RAMPS
SKU: 99AAC1600100



CITY 2 PLUS - CUSHION SETS
GREEN SKU: 22ACM0100103
PURPLE SKU: 22ACM0100105
BLUE SKU: 22ACM0100102
RED SKU: 22ACM0100106