Automatic Blood Pressure Monitor, Wrist Model



HCPCS: A4670-1,A4670-1



FIGURE A

- Large easy-to-read digital display
- 120 memory recall
- Date & time display
- Pulse Rate Indicator
- Auto-off function
- Lightweight & easy to use
- Quiet inflation method
- Wrist model

PRODUCTS	Description
BP2116	BP Monitor Wrist Auto 1 EA

SPECIFICATIONS	Description
Warranty	2 Years
Measurement Range	Pressure: 0~300mmHg +/-3mmHg
Power Supply	2 x AAA Batteries