
P703



1. Hi/Lo Motor- SP01-TA23-1060-001
   Hi/Lo & Foot Motor Extension Cord- SP01-P703-80E
2. Head Motor- SP01-TA31-1060-001
3. Foot Motor- SP01-TA31-1060-002
4. Motor Control Box- SP01-TC21-1060-001
   "Y" Cable- SP01-P703-70E
   Power Cord- SP01-P703-60-02E
5. Hand Pendant (8 button)- SP02-TH21060003
6. Head Deck Truck Assembly Left- SP01-P703-04-03E
7. Head Deck Truck Assembly Right- SP01-P703-04-02E
8. Foot Deck Truck Assembly- SP01-P703-20E
9. Head Deck w/Hardware- SP01-P703-10001E
10. Head Board Adjustment Bracket- SP01-P703-10008E
11. Head Board Mounting Bracket- SP01-PCB901-10003E
12. Head/Foot Board Mounting Bolt (ea)- SP01-UF16108061E
13. Head End Mattress Retainer (ea)- SP01-P703-02-03E
14. Width Mattress Retainer (ea)- SP01-P703-02-02E
15. Head/Thigh Deck Hardware- SP01-P703-10002E
16. Thigh Deck w/Hardware- SP01-P703-10003E
17. Foot Deck w/Hardware- SP01-P703-10005E [3]
    Foot Deck w/Hardware- SP01-P703-10005AE [4]
    Foot Deck Hardware- SP01-P703-10004E
18. Foot Board Adj. Length Mtg. Bracket w/Hardware- SP01-P703-10009E [3]
    Foot Board Adj. Length Mtg. Bracket w/Hardware- SP01-P703-10009AE [4]
19. Foot Board Mounting Bracket w/Hardware (pr)- SP01-PCB901-10019E
20. Foot End Mattress Retainer- SP01-P703-03-03E
21. Foot Deck Ratcheting Assembly- SP01-PCB300-01-57E
22. Wall Guard- SP01-P703-30E
23. Brake Stop- SP01-PCB900-24BS
24. Locking Caster- SP01-TL500-37E
25. Non-Locking Caster- SP01-TL500-29E
26. End Cap- PC20303001E-CO

Not Shown

Head Width Extension Left- SP01-P703WE-10002E
Head Width Extension Left- SP01-P703WE-2EAI [1]
Head Width Extension Right- SP01-P703WE-10001E
Head Width Extension Right- SP01-P703WE-1EAI [2]
Thigh Width Extension Left- SP01-P703WE-10004E
Thigh Width Extension Right- SP01-P703WE-10003E
Foot Width Extension Left- SP01-P703WE-10006E
Foot Width Extension Right- SP01-P703WE-10005E
Bed Angle Indicator Left- SP01-ANGLE-L
Bed Angle Indicator Right- SP01-ANGLE-R

[1] For use SP01-ANGLE-L
[2] For use SP01-ANGLE-R
[3] 5/2023 and before.
[4] 6/2023 and later.

12/15/2025