ZooMe Auto Flex



1. Controller- FLEX-29
2. Throttle Pot- FLEX-AUTO-01
2. Throttle Pot Assembly- FLEX-AUTO-02
   Throttle Lever- FLEX-AUTO-03
3. On-Off Switch- FLEX-AUTO-05
4. Speed Pot- FLEX-AUTO-06
5. Horn Button- FLEX-AUTO-07
   Horn- FLEX-AUTO-04
6. Junction Board- FLEX-AUTO-08 [7]
   Junction Board- FLEX-AUTO-08N [2]
7. Battery Gauge- FLEX-AUTO-09
8. Upper Console Blue- FLEX-AUTO-10
   Upper Console Red- FLEX-AUTO-10R
9. Lower Console Blue- FLEX-AUTO-11
   Lower Console Red- FLEX-AUTO-11R
10. Quick Release Clamp- FLEX-AUTO-12
11. Handle Connector- FLEX-AUTO-13
12. Steering Knuckle, Left- FLEX-AUTO-14
    Steering Knuckle, Right- FLEX-AUTO-15
    Steering Knuckle Bearing- FLEX-11
13. Footboard Blue- FLEX-AUTO-16
    Footboard Red- FLEX-AUTO-16R
14. Floor Mat- FLEX-AUTO-17
15. Folding Actuator Assembly- FLEX-AUTO-18
    Folding Actuator Assy- FLEX-AUTO-51 [1]
    Folding Actuator Assy- FLEX-AUTO-51-2.4 [5]
16. Wireless Receiver- FLEX-AUTO-19
    Wireless Receiver- FLEX-AUTO-50 [1]
    Wireless Receiver- FLEX-AUTO-50N [2]
    Wireless Receiver- FLEX-AUTO-50-2.4 [5]
17. Key- FLEX-AUTO-20
    Wireless Remote Key- FLEX-AUTO-21
    Wireless Remote Key- FLEX-AUTO-21-2.4 [5]
18. Armrest Assembly, Left- FLEX-AUTO-23
    Armrest Assembly, Right- FLEX-AUTO-24
    Arm Pad- FLEX-37
19. Armrest Knob- FLEX-AUTO-25
20. Charger Port- FLEX-AUTO-26
21. Actuator Switch- FLEX-AUTO-27
    Actuator Switch- FLEX-AUTO-27N [8]
    Actuator Switch- FLEX-AUTO-27-2.4G [5]
22. Reverse Horn- FLEX-AUTO-28
23. Rear Shroud Blue- FLEX-AUTO-29
    Rear Shroud Red- FLEX-AUTO-29R
    Rear Shroud Decor Blue- FLEX-AUTO-30
    Rear Shroud Decor Red- FLEX-AUTO-30R
24. Anti-Tipper Wheel- FLEX-AUTO-31
25. Main Harness- FLEX-AUTO-45
    Main Harness- FLEX-AUTO-45N [2]
    Main Harness- FLEX-AUTO-45-2.4 [5]
26. Front Shroud Blue- FLEX-12
    Front Shroud Red- FLEX-12R
27. Gas Spring Assembly- FLEX-AUTO-34
28. Front Reflector Assy Blue- FLEX-AUTO-39
    Front Reflector Assy Red- FLEX-AUTO-39R
29. Front Wheel- FLEX-07
30. Release Cable- FLEX-AUTO-46
31. Release Lever- FLEX-21
32. Battery Pack- [3] [6]
    Battery Pack- FLEX-AUTO-47 [4]
    Battery Case Only- FLEX-AUTO-48 [4] [7]
    Battery Contacts (4 ea)- FLEX-31
33. Rear Wheel- FLEX-25
34. Battery Charger- FLEX-27N
35. Motor/Brake/Transaxle Assy- FLEX-AUTO-44
    Transaxle- FLEX-51
    Motor- FLEX-52
    Brake- FLEX-AUTO-43
36. Seat Complete- FLEX-AUTO-37
    Back Upholstery- FLEX-AUTO-36
    Seat Upholstery- FLEX-AUTO-35
37. Tie Rod Complete- FLEX-13
38. Hub Cap- FLEX-10
39. Wheel/Stem Bearing- FLEX-11
40. Armrest Receiver, Left- FLEX-AUTO-32
    Armrest Receiver, Right- FLEX-AUTO-33
41. Folding Latch Assy Yellow- FLEX-AUTO-38
    Folding Latch Assy Red- FLEX-AUTO-49 [1]

not shown Ext Charging Station- FLEX-50
not shown Basket Only- FLEX-AUTO-40
not shown Basket mount- FLEX-AUTO-41
not shown Cup Holder- FLEX-AUTO-42

[1] Serial # (21)24S1803000001 to (21)24S1909000001
[2] Serial # (21) 24S1909000001 and higher
[3] Lower than Serial # (21)24S1712000001
[4] Serial # (21)24S1712000001 & higher
[5] Serial #(21)24S2112000001 & higher
[6] No longer available.
[7] When stock is gone the part is no longer available.
[8] Serial #24S1803000001 to (21)24S2112000001

12/02/2024