Scout Compact 3-Wheel Scooter



HCPCS: K0800-1,K0800


FIGURE A


FIGURE B


FIGURE C


FIGURE D


FIGURE E


FIGURE F

The Scout transportable scooter combines easy and convenient mobility with superior value. The Scout is compact, lightweight and easy to disassemble-perfect for taking wherever you want to go, either around the neighborhood or far from home. It can be disassembled quickly and easily into five easy-to-handle pieces for loading into a car, truck, or van and taking shopping or on vacation. The three-wheel version of the Scout allows for nimbly negotiating indoor obstacles. The Scout also offers great features like a quick connect battery pack and a delta-style tiller to make operating the scooter easier for those with limited strength or dexterity. This model comes with color-through panels to hide scuffs and scratches and black wheel rims.

- Quick and easy disassembly
- Interchangeable color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy freewheel operation

PRODUCTS	Description
SFSCOUT3	Scout 3-Wheel Scooter, Red/Blue, 1/ea

SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	85 lb
Product Dimensions	42.5" (L) x 20.5" (W) x 34" (H)
Seat Dimensions	16" (D) x 16.5" (W)
Weight Capacity	300 lb
Top Speed	4.25 mph
Batteries (pair)	12V x 12AH
Battery Charger	2A Offboard
Climbing Angle	6°
Color	Red & Blue
Controller	Dynamic 45A
Ground Clearance	2.5"

SPECIFICATIONS	Description
HCPCS	K0800
Maximum Range	9 miles
Motor	24V x 270W x 4700 RPM
Seat Type	Stadium
Turning Radius	45.5"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 8" x 2" | Rear: 8" x 2"